

User Experience Upgrade

General Changes

- Navigation menu is now collapsible.
- New color-coded notification cards have been added to display feedback information.
- Increased the number of athletes shown in athlete selection windows.
- Support for tooltips has been added to give users relevant information.
- Notification when application loses and regain internet connection.

Athlete Screen

- Redesigned athlete basic info page.
- Athlete number is no longer required to be unique.
- The new profile image selector now lets the user crop the image in situ without redirecting to an additional modal.
- The measurement page now shows all inputs' minimum and maximum value ranges.
- The Dominance page is now color-coded to make the chosen side more straightforward.
- Athlete Dashboard charts now support more touch interactions. Such as pinch, zoom and swipe to move chart view.
- New chart controls for charts in the Athlete Dashboard have been added.
- Deep KT Chart, Superficial KT Chart, and Strength KT Chart have been retired.
- Athlete Region Overview cards are now color-coded to match injury alerts.
- Users can now move between athletes by pressing the athlete profile image and using the popup modal.

Planning Screen

- Move and Copy functionality have been made visible.
- Move and Copy buttons in the parent card can

be used to move/copy sessions from one week to another or from one month to another.

Testing & Training Screens

- Users can now select any session that shares the same selected athlete and date.
- The protocol description modal layout has been improved.
- Protocol Image dimensions have now been constrained for tablet viewing.
- Protocol Setup and Cue text contrast have been increased to improve readability.
- The protocol angle and unit height modal layout have been improved.
- Protocol angle controls now give the user more feedback on target angles.
- Inclusion of a loading animation to force trace when the side is switched to hide interference from KT360.
- Inclusion of a tooltip for connected sensor information.
- Increased number of warnings for unsaved data.
- The x-axis of the post-session chart has been corrected.
- In training, can no longer delete reps before the current set.
- Text filter added to template list so that templates can be found quickly.
- Improved functionality for Calf Capacity test (eg. BPM adjustment, manual input of reps).

Testing & Training Session Listing

- "Forward" and "Back" buttons were added to the date picker to increment through dates one day at a time.
- Pressing "Resume Session" in the planning screen will open the selected session in the test session listing screen.

- The athlete search function will no longer show empty sessions if the chosen athlete is not found in said sessions.
- Can now move between dates with an athlete entered in the search function and have the sessions filtered based on the selected athlete.

Reports

- Added Pre-selected Parameters for all reports.
- In-report regeneration functionality has been improved.
- Users can now edit all parameters from within the generated report. The user no longer needs to exit the report to change protocol.
- Group Rank report normalization can now be changed from within a generated report.
- Users can now see a list of athletes with valid monitoring data from the generation modal.
- Increased contrast to improve report readability for all reports.

Account Screen

- No longer display user email in the user list.
- Added a user preference to turn off all time max target line.
- Added Color-coded buttons depending on function.
- Added warnings before removing admin privileges.
- Added warning before removing Multi-Factor Authentication.

Export Screen

- Added new error messages for failed exports.
- Added new layout for future improvements.